

Product Datasheet

Mouse GPX4 ELISA Kit (orb1479370)

Description	Mouse Glutathione Peroxidase 4 (GPX4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids; sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	070325
Sensitivity	0.056ng/mL
Entrez	625249
Expiration Date	Please enquire.