

Product Datasheet

Human CLDN2 ELISA Kit (orb1479390)

Description	Human Claudin 2 (CLDN2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: tissue homogenates, cell lysates and other biological fluids; sandwich
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P57739
Sensitivity	0.052ng/mL
Entrez	9075
Expiration Date	Please enquire.