

Product Datasheet

Progesterone Receptor (Phospho-Ser190) Antibody (orb14799)

Description Progesterone Receptor (Phospho-Ser190)
Antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note For research use only


Clonality Recombinant

Uniprot ID **P06401**


Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of extract from SK...



Immunofluorescence staining of methanol-...