

Product Datasheet

Tsp23 antibody [Out of stock] (orb147992)

Description

Mouse monoclonal to Tsp23 antibody

Species/Host	Mouse
Conjugation	FITC
Tested Applications	AM, IHC, WB
Immunogen	Recombinant tsp23
Target	TSP23
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-194 was sufficient for detection of tsp23 in 20 µg of transiently transfected Hela cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1 Kappa
Clonality	Monoclonal
MW	19kDa
Uniprot ID	B9A003
NCBI	NP_001136126.1
Entrez	80755
Dilution Range	WB (1:1000), IHC (1:1000)
Expiration Date	12 months from date of receipt.