



Product Datasheet

PKD/PKC μ (Phospho-Ser738) Antibody (orb14803)

Description	PKD/PKC μ (Phospho-Ser738) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q15139
Expiration Date	12 months from date of receipt.



Western blot analysis of extract from A4...



Immunofluorescence staining of methanol...