

Product Datasheet

MDC1 antibody [Out of stock] (orb148035)

Description

Mouse monoclonal to MDC1 (ATTO-594).

checkpoint protein 1, plays a role...

Species/Host

Mouse

Reactivity

Bovine, Human, Mouse, Primate

Conjugation

ATTO-594

Tested

ICC, WB

Applications

Immunogen

GST-tagged recombinant protein corresponding to mouse MDC1 at and around the N-terminus

Target

MDC1

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.5 µg/ml of SMC-197 was sufficient for detection of MDC1 in 10 µg of HeLa cell lysate by ECL immunoblot analysis.

Isotype

IgG1

Clonality

Monoclonal

MW

184kDa

Uniprot ID

[Q5PSV9](#)

NCBI

[NP_001010833.2](#)

Entrez

[240087](#)

Dilution Range

WB (1:2000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.