

Product Datasheet

EndoPDI antibody [Out of stock] (orb148135)

Description

Mouse monoclonal to EndoPDI (ATTO-

isomerase (EndoPDI) is a thior...

Species/Host

Mouse

Reactivity

Human

Conjugation

ATTO-488

Tested

ELISA, ICC, IHC, IP, WB

Applications

Immunogen

Synthesized peptide - 12 amino acids long (ADGEDGQDPHSK) corresponding to residues 52 - 63 of the EndoPDI protein

Target

PDI

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-204 was sufficient for detection of EndoPDI in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti - mouse IgG:HRP as the secondary antibody.

Isotype

IgG2a

Clonality

Monoclonal

MW

48kDa

Uniprot ID

[Q8NBS9](#)

NCBI

[NP_001139021](#)

Entrez

[81567](#)

Dilution Range

WB (1:200)

Expiration Date

12 months from date of receipt.