



Product Datasheet

p53 (Phospho-Ser6) Antibody (orb14815)

Description	p53 (Phospho-Ser6) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P04637
Expiration Date	12 months from date of receipt.



Western blot analysis of p53(Phospho-Ser...



Immunohistochemical analysis of formalin...