


Product Datasheet

Androgen Receptor (Phospho-Ser650) Antibody (orb14842)

Description	Androgen Receptor (Phospho-Ser650) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P10275
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...