

Product Datasheet

TrpV3 antibody (orb148488)

Description

Mouse monoclonal to TrpV3 (PerCP). Ion

	help establish a...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	PerCP
Tested Applications	AM, ICC, IHC, IP, WB
Immunogen	Synthetic peptide amino acids 774-791 (C-terminus) of rat TrpV3
Target	TRPV3
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-319 was sufficient for detection of TrpV3 in 10 µg of COS-1 cell lysate transiently transfected with TrpV3 by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2a
Clonality	Monoclonal
MW	70kDa
Uniprot ID	Q4QYD9
NCBI	NP_001020928
Entrez	497948
Dilution Range	WB (1:1000), IHC (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.