



Product Datasheet

NMDAR2B (phospho-Tyr1474) Antibody (orb14881)

Description	NMDAR2B (phospho-Tyr1474) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13224
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of methanol-...



Western blot analysis of NMDAR2B (phospho-...