

## Product Datasheet

### GABA-A Receptor Delta antibody [Out of stock] (orb148885)

**Description**

Mouse monoclonal to GABA-A R Delta

a member of the superfamily of...

**Species/Host**

Mouse

**Reactivity**

Human, Mouse, Rat

**Conjugation**

ATTO-594

**Tested**

AM, ICC, WB

**Applications**
**Immunogen**

Synthetic peptide amino acids 15-34 (N-terminus) of rat GABA-A-R-Delta

**Target**

GABA A Receptor

**Preservatives**

PBS pH 7.4, 50% glycerol, 0.09% sodium azide \*Storage buffer may change when conjugated

**Concentration**

1 mg/ml

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

2 µg/ml of SMC-345 was sufficient for detection of Delta1 GABA-A receptor in 10 µg of rat brain lysate by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.

**Isotype**

IgG2a

**Clonality**

Monoclonal

**MW**

55kDa

**Uniprot ID**

[P18506](#)

**NCBI**

[NP\\_058985.1](#)

**Entrez**

[29689](#)

**Dilution Range**

WB (1:1000), IHC (1:1000), ICC/IF (1:100)

**Expiration Date**

12 months from date of receipt.