



Product Datasheet

ADD1 (Phospho-Ser726) Antibody (orb14892)


Description	ADD1 (Phospho-Ser726) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P35611
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of ADD1(Phospho-Se...



Confocal immunofluorescent analysis of H...