


Product Datasheet

Paxillin (Phospho-Tyr31) Antibody (orb14904)

Description	Paxillin (Phospho-Tyr31) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P49023
Expiration Date	12 months from date of receipt.



Western blot analysis of Paxillin(Phospho-Tyr31)