

Product Datasheet

Aldehyde Dehydrogenase Microplate Assay Kit [Out of stock] (orb1494272)

Description	Detection and Quantification of Aldehyde Dehydrogenase Activity.
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
CAS Number	9028-86-8
Expiration Date	12 months from date of receipt.