



## Product Datasheet

# p70 S6 Kinase (Phospho-Thr421) Antibody (orb14943)


|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | p70 S6 Kinase (Phospho-Thr421) Antibody                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                    |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | IF, IHC, WB                                                                               |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                     |
| <b>Clonality</b>           | Recombinant                                                                               |
| <b>Uniprot ID</b>          | <b>P23443</b>                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |



Western blot analysis of p70 S6 Kinase(P...



Immunohistochemical analysis of formalin...



Confocal immunofluorescent analysis of M...