

Product Datasheet

Human IL-21 ELISA Kit [Out of stock] (orb1494460)

Description	Sandwich ELKSA for quantitative detection of Human IL-21 concentrations.
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL21
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9HBE4
Sensitivity	15 pg/ml
Entrez	59067
Expiration Date	Please enquire.