

Product Datasheet

Human IL-1RL1 ELISA Kit [Out of stock] (orb1494461)

Description	Sandwich ELKSA for quantitative detection of Human IL-1RL1 concentrations.
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL1RL1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q01638
Sensitivity	15 pg/ml
Entrez	9173
Expiration Date	Please enquire.