

Product Datasheet

Human IL-18 BP α ELISA Kit [Out of stock] (orb1494463)

Description	Sandwich ELKSA for quantitative detection of Human IL-18 BP α concentrations.
Range	93.6 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	IL18BP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell culture supernates, Serum, Plasma
Uniprot ID	095998
Sensitivity	50 pg/ml
Entrez	10068
Expiration Date	Please enquire.