

Product Datasheet

Human IL-17E ELISA Kit [Out of stock] (orb1494464)

Description	Sandwich ELKSA for quantitative detection of Human IL-17E concentrations.
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	IL25
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9H293
Sensitivity	10 pg/ml
Entrez	64806
Expiration Date	Please enquire.