

Product Datasheet

Human IL-12 p70 ELISA Kit [Out of stock] (orb1494467)

Description	Sandwich ELKSA for quantitative detection of Human IL-12 p70 concentrations.
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL12A; IL12B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P29460, P29459
Sensitivity	7 pg/ml
Entrez	3592, 3593
Expiration Date	Please enquire.