

Product Datasheet

PIEZO1 polyclonal antibody (orb1495188)

Description	PIEZO1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length protein of human PIEZO1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 286 kDa
Uniprot ID	Q92508
Dilution Range	WB (1/500 - 1/2000)
Expiration Date	12 months from date of receipt.