

Product Datasheet

Thyroglobulin monoclonal antibody (orb1495288)

Description	Thyroglobulin monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Thyroglobulin. The exact sequence is proprietary.
Preservatives	Mouse IgG1. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 304 kDa
Uniprot ID	P01266
Dilution Range	IHC (1/100 - 1/300)
Expiration Date	12 months from date of receipt.