

Product Datasheet

Beta2A-tubulin monoclonal antibody (orb1495616)

Description	Beta2A-tubulin monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Beta2A-tubulin. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 50 kDa
Uniprot ID	Q13885
Dilution Range	WB (1/50000 - 1/100000), IP (1/100 - 1/200)
Expiration Date	12 months from date of receipt.