

Product Datasheet

Acetyl Lysine monoclonal antibody (orb1495629)

Description	Acetyl Lysine monoclonal antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Purified protein corresponding to Acetyl Lysine.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB (1/1000 - 1/2000), IHC (1/200 - 1/500)
Expiration Date	12 months from date of receipt.