

Product Datasheet

CD45RB Rat monoclonal antibody (orb1495666)

Description	CD45RB Rat monoclonal antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Cloned mouse Th2 cell lines
Preservatives	Rat IgG2a kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 144 kDa
Uniprot ID	P06800
Dilution Range	IF (1/50 - 1/200)
Expiration Date	12 months from date of receipt.