

Product Datasheet

Kir5.1 (phospho-S416) polyclonal antibody (orb1495829)

Description	Kir5.1 (phospho-S416) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Kir5.1 with a phosphorylation site at S416. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	Q9NPI9
Dilution Range	WB (1/500 - 1/1000), IHC (1/50 - 1/200), IF/ICC (1/50 - 1/200)
Expiration Date	12 months from date of receipt.