

Product Datasheet

IκB-β (Phospho-S23) polyclonal antibody [KU] (orb1496028)

Description	IκB-β (Phospho-S23) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human IκB-β around the phosphorylation site of Serine 23.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 42 kDa
Uniprot ID	Q15653
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.