

Product Datasheet

ZO-1 polyclonal antibody (orb1496077)

Description	ZO-1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-273 of human ZO-1 (NP_003248.3).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 250 kDa
Uniprot ID	Q07157
Dilution Range	WB: 1:500~1:2000 IF: 1:50~1:200
Expiration Date	12 months from date of receipt.