



Product Datasheet

PKR (Phospho-Thr446) Antibody (orb14961)


Description	PKR (Phospho-Thr446) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P19525
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of PKR(Phospho-Thr...



Confocal immunofluorescent analysis of H...