

Product Datasheet

Mouse uPAR / PLAUR Protein (orb1496133)

Description	Urokinase plasminogen activator surface receptor (U-PAR) is also known as PLAUR, Monocyte activation...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	uPAR
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 31.7 kDa. The protein migrates as 45-63 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_035243.1
Purity	85%
Source	Mouse uPAR, His Tag (orb1496133) is expressed from human 293 cells (HEK293). It contains AA Leu 24 - Pro 296 (Accession # P35456-1).
MW	31.7 kDa
NCBI	NP_035243.1
Expiration Date	6 months from date of receipt.