

Product Datasheet

Biotinylated Human Complement C5 Protein (orb1496162)

Description

Biotinylated Human Complement C5, His, is

Biotinylated Human Complement C5, His, is
0.2 µg of Biotinylated Human Complement C5, His, per well

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

Complement C5

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;
-70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 189.9 kDa. As a result of glycosylation, the protein migrates as 75 kDa, 120 kDa and 200 kDa under reducing (R) condition, and 200 kDa under non-reducing (NR) condition (SDS-PAGE).

Purity

95%

Source

Biotinylated Human Complement C5, His, (orb1496162) is expressed from human 293 cells (HEK293). It contains AA Gln 19 - Cys 1676 (Accession # P01031-1).

MW

189.9 kDa

Uniprot ID


P01031

NCBI

NP_001726.2

Expiration Date

6 months from date of receipt.



Immobilized
Biotinylated
Human
Complement C5