

Product Datasheet

Mouse IL-1RL1 / ST2 Protein (orb1496167)

Description

Mouse IL-1RL1, Fc Tag is expressed from

Ser 27 - Arg 332 (A...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-1RL1

Preservatives

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 61.5 kDa. The protein migrates as 80-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Mouse IL-1RL1, Fc Tag (orb1496167) is expressed from human 293 cells (HEK293). It contains AA Ser 27 - Arg 332 (Accession # P14719-1).

MW

61.5 kDa

Uniprot ID

P14719

Expiration Date

6 months from date of receipt.