


Product Datasheet

Human IL-10 Protein (orb1496180)

Description	Human IL-10, premium grade is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Asn 1...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	IL-10
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 18.6 kDa. The protein migrates as 17 kDa±3 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-10, premium grade (orb1496180) is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Asn 178 (Accession # P22301-1).
MW	18.6 kDa
NCBI	NP_000563.1
Expiration Date	6 months from date of receipt.



Immobilized
Human IL-10,
premium
grade (...)