

Product Datasheet

Human CCL24/Eotaxin-2 Protein (orb1496212)

Description

Human CCL24, His Tag is expressed from

Val 27 - Cys 119 (Ac...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CCL24

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 12.3 kDa. The protein migrates as 18-21 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human CCL24, His Tag (orb1496212) is expressed from human 293 cells (HEK293). It contains AA Val 27 - Cys 119 (Accession # O00175).

MW

12.3 kDa

Uniprot ID
O00175
Expiration Date

6 months from date of receipt.