


Product Datasheet

Cynomolgus IL-10 Protein, His Tag (orb1496244)

Description	Interleukin-10 (IL-10) is also known as human cytokine synthesis inhibitory factor (CSIF), is an ant...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-10
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 20.6 kDa. The protein migrates as 21-22 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Cynomolgus IL-10 Protein, His Tag (orb1496244) is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Asn 178 (Accession # P79338).
MW	20.6 kDa
Expiration Date	6 months from date of receipt.



Immobilized
Cynomolgus
IL-10 R
alpha Pro...