

Product Datasheet

Human Oncostatin M (OSM) (26-221) Protein, Tag Free (MALS verified) (orb1496309)

Description

Oncostatin M is also known as

belonging to the interleukin-6 family of cytoki...

Reactivity

Human

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

Oncostatin M

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 22.2 kDa. The protein migrates as 25-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human Oncostatin M (OSM) (26-221), Tag Free (orb1496309) is expressed from human 293 cells (HEK293). It contains AA Ala 26 - Arg 221 (Accession # P13725).

MW

22.2 kDa

Expiration Date

6 months from date of receipt.