


Product Datasheet

Human IL-10 R alpha / CD210 Protein (orb1496310)

Description

Human IL-10 R alpha Protein, Fc Tag is expressed

	2...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-10 R alpha
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 50.8 kDa. The protein migrates as 65-75 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-10 R alpha Protein, Fc Tag (orb1496310) is expressed from human 293 cells (HEK293). It contains AA His 22 - Asn 235 (Accession # Q13651).
MW	50.8 kDa
Uniprot ID	Q13651
Expiration Date	6 months from date of receipt.



Loaded
Human IL-
10 R alpha /
CD210, Fc
T...