



Product Datasheet

p70 S6 Kinase (Phospho-Ser424) Antibody (orb14964)


Description	p70 S6 Kinase (Phospho-Ser424) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P23443
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of p70 S6 Kinase(P...



Confocal immunofluorescent analysis of H...