

## Product Datasheet

### Zebrafish IL8 ELISA Kit (orb1496456)

|                        |                                     |
|------------------------|-------------------------------------|
| <b>Description</b>     | Zebrafish IL8 ELISA Kit             |
| <b>Species/Host</b>    | Zebrafish                           |
| <b>Reactivity</b>      | Zebrafish                           |
| <b>Range</b>           | 7.82-500 pg/mL                      |
| <b>Note</b>            | For research use only               |
| <b>Sample Types</b>    | serum, plasma and biological fluids |
| <b>Sensitivity</b>     | 3.1 pg/mL                           |
| <b>Expiration Date</b> | Please enquire.                     |

