

Product Datasheet

Rat NACHT, LRR and PYD domains-containing protein 6 ELISA Kit (orb1496496)

Description	Rat NACHT, LRR and PYD domains-containing protein 6 ELISA Kit
Species/Host	Rat
Reactivity	Rat
Range	0.79-50 ng/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids
Uniprot ID	Q63035
Sensitivity	0.28 ng/mL
Expiration Date	Please enquire.

