

Product Datasheet

VGLUT1 antibody (orb149726)

Description

Mouse monoclonal to VGLut1 (FITC). VGLUT1 is

and transpor...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested

IHC, WB

Applications

Immunogen

Fusion protein amino acids 493-560
(cytoplasmic C-terminus) of rat VGLut1

Target

VGLUT1

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-394 was sufficient for detection
of VGLut1 in 20 µg of rat brain lysate by
colorimetric immunoblot analysis using goat
anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

52kDa

Uniprot ID

[Q62634](#)

NCBI

[NP_446311.1](#)

Entrez

[116638](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.