


Product Datasheet

Keratin 8 (Phospho-Ser74) Antibody (orb14976)

Description	Keratin 8 (Phospho-Ser74) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05787
Expiration Date	12 months from date of receipt.



Western blot analysis of Keratin 8(Phosp...