

Product Datasheet

VGLUT3 antibody [Out of stock] (orb149769)

Description

Mouse monoclonal to VGlut3 (ATTO-594).

(VGLUT3) is a multipass mem...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-594

Tested

ICC, IHC, WB

Applications
Immunogen

Fusion protein amino acids 546-588
(cytoplasmic C-terminus) of rat VGLUT3

Target

VGLUT3

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium
azide *Storage buffer may change when
conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored
at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-397 was sufficient for
detection of VGLUT3 in 20 µg of CV-1
fibroblast cells (lysate) transfected with
VGLUT3 by colorimetric immunoblot
analysis using goat anti-mouse IgG:HRP
as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

65kDa

Uniprot ID

[Q7TSF2](#)

NCBI

[NP_714947.1](#)

Entrez

[266767](#)

Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.