



Product Datasheet

Shc1 (Phospho-Tyr349) Antibody (orb14980)

Description	Shc1 (Phospho-Tyr349) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P29353
Expiration Date	12 months from date of receipt.



Western blot analysis of Shc1(Phospho-Ty...



Confocal immunofluorescent analysis of H...