



Product Datasheet

Shc1 (Phospho-Tyr427) Antibody (orb14981)


Description	Shc1 (Phospho-Tyr427) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P29353
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescent
analysis of H...



Immunohistochemical
analysis of formalin...



Western blot analysis
of Shc1(Phospho-Ty...