

Product Datasheet

Ankyrin B antibody (orb149831)

Description

Mouse monoclonal to Ankyrin B (PerCP).

mediate the at...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

PerCP

Tested

ICC, IHC, WB

Applications

Immunogen

Synthetic peptide amino acids 203-496 of human Ankyrin-B

Target

Ankyrin B

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-400 was sufficient for detection of Ankyrin-B in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

200kDa

Uniprot ID

[Q01484](#)

NCBI

[NP_001120965.1](#)

Entrez

[287](#)

Dilution Range

WB (1:1000), IHC (1:200), ICC/IF (1:200)

Expiration Date

12 months from date of receipt.