

Product Datasheet

SNAT1 antibody [Out of stock] (orb149835)

Description

Mouse monoclonal to SNAT1 (ATTO-488).

	transporters (SNAT) of t...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-488
Tested Applications	FC, ICC, IHC, WB
Immunogen	Fusion protein amino acids 1- 63 of rat SNAT1
Target	SLC38A1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-401 was sufficient for detection of SNAT1 in 20 µg of lysates from neocortical neurons cultured under amino acid starvation conditions and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	50kDa
Uniprot ID	Q9JM15
NCBI	NP_620187.1
Entrez	170567
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.