

Product Datasheet

GABA-B Receptor 1 antibody [Out of stock] (orb149871)

Description

Mouse monoclonal to GABA-B R 1

aminobutyric acid) is the primary inhibitory ...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-594

Tested

ICC, WB

Applications
Immunogen

Fusion protein amino acids 873-977 (cytoplasmic C-terminus) of rat GABA(B)R1

Target

GABA B Receptor 1

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-403 was sufficient for detection of GABA(B)R1 in 20 µg of rat brain membrane lysate and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

115kDa

Uniprot ID

[Q9Z0U4](#)

NCBI

[NP_112290.2](#)

Entrez

[81657](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.