

Product Datasheet

GABA-B Receptor 1 antibody (orb149878)

Description

Mouse monoclonal to GABA-B R 1 (Biotin).

	primary inhibitory ne...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, WB
Immunogen	Fusion protein amino acids 873-977 (cytoplasmic C-terminus) of rat GABA(B)R1
Target	GABA B Receptor 1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-403 was sufficient for detection of GABA(B)R1 in 20 µg of rat brain membrane lysate and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	115kDa
Uniprot ID	Q9Z0U4
NCBI	NP_112290.2
Entrez	81657
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.